

Rabbit Anti-Human E-Cadherin/Cd324 Monoclonal Antibody For Flow Cytometry - P12830

Item Number: CM0009559



Introduction

-conjugated rabbit monoclonal antibody targeting human E-cadherin/CD324, validated for flow cytometry. Ideal for studying epithelial cell adhesion, EMT, cancer invasion, and cadherin-mediated cell sorting. Reactivity: human; molecular weight: ~90–97 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human E-Cadherin/CD324 mAb
Remarks/Alias	UVO; CDHE; ECAD; LCAM; Arc-1; BCDS1; CD324
Species	Human
Gene ID	999
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 155-707 of human E-Cadherin/CD324 (NP_004351.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Human
SWISS	P12830
Protein Weight	90kDa/97kDa
Shipping	Ice bag