

Rabbit Anti-Human Egfr Mab For Fc - P00533

Item Number: CM0008412



Introduction

Rabbit anti-Human EGFR monoclonal antibody validated for flow cytometry in human samples. Targets extracellular domain (aa 335-525), enabling receptor signaling studies and oncology research applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human EGFR mAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2
Species	Human
GeneID (Human)	1956
Immunogen	Recombinant protein\Recombinant fusion protein containing a sequence corresponding to amino acids 335-525 of human EGFR (NP_005219.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P00533
Protein Weight	44kDa/69kDa/77kDa/134kDa
Shipping	Ice bag