

Rabbit Anti-Human Egfr Monoclonal Antibody For If/Icc, Fc - P00533

Item Number: CM0004740



Introduction

-conjugated rabbit anti-human EGFR monoclonal antibody for immunofluorescence and flow cytometry. Targets the extracellular domain (aa 335-525) of the epidermal growth factor receptor. Suitable for human samples; validated for research use.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human EGFR mAb
Remarks/Alias	EGFR; ERBB; ERBB1; HER1; NISBD2; PIG61; mENA; epidermal growth factor receptor
Species	Human
GeneID (Human)	1956
GeneID	1956
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 335-525 of human EGFR (NP_005219.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	P00533
Protein Weight	44kDa/69kDa/77kDa/134kDa
Shipping	Ice bag