


Rabbit Anti-Human Erbb3/Her3 Monoclonal Antibody () For Flow Cytometry - P21860

Item Number: CM0011109



Introduction

Rabbit monoclonal antibody conjugated to targeting human ErbB3/HER3 (UniProt P21860). Validated for flow cytometry (FC). Ideal for studying ErbB3/HER3 expression and neuregulin signaling in human epithelial and cancer samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ErbB3/HER3 mAb
Remarks/Alias	ERBB3; ErbB-3; HER3; LCCS2; MDA-BF-1; c-erbB-3; c-erbB3; erbB3-S; p180-ErbB3; p45-sErbB3; p85-sErbB3; erb-b2 receptor tyrosine kinase 3
Species	Human
Gene ID (Human)	2065
Gene ID	2065
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-643 of human ErbB3/HER3 (NP_001973.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21860
Protein Weight	20kDa/36kDa/77kDa/141kDa/148kDa
Shipping	Ice bag