

Rabbit Anti-Human Erbb3/Her3 Monoclonal Antibody For Flow Cytometry - P21860

Item Number: CM0010819



Introduction

-conjugated rabbit monoclonal antibody targeting human ErbB3/HER3 for flow cytometry. Recognizes the extracellular domain (aa20-643) of the receptor tyrosine kinase involved in neuregulin signaling and PI3K activation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ErbB3/HER3 mAb
Remarks/Alias	HER3; FERLK; LCCS2; VSCN1; ErbB-3; c-erbB3; erbB3-S; MDA-BF-1; c-erbB-3; p180-ErbB3; p45-sErbB3; p85-sErbB3
Species	Human
Gene ID	2065 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-643 of human ErbB3/HER3 (NP_001973.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (FC)
Cross Reactivity	Human
Swiss-Prot	P21860
Protein Weight	148 kDa
Shipping	Ice bag