

Rabbit Anti-Human Erg Monoclonal Antibody For Flow Cytometry (Intracellular) - P11308

Item Number: CM0010590



Introduction

-conjugated rabbit monoclonal antibody targeting human ERG (p55/erg-3). Specific for intracellular flow cytometry applications, validated in human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ERG mAb
Remarks/Alias	p55; erg-3
Species	Human
GeneID (Human)	2078
GeneID	2078
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human ERG (NP_891548.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P11308
Protein Weight	54kDa
Shipping	Ice bag