

Rabbit Anti-Human Erg Monoclonal Antibody For Intracellular Flow Cytometry - P11308

Item Number: CM0008692



Introduction

-conjugated rabbit monoclonal antibody targeting human ERG (p55/erg-3). Designed for intracellular flow cytometry (FC intra) applications. Recognizes human ERG; predicted molecular weight 54kDa. UniProt: P11308.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ERG mAb
Remarks/Alias	p55; erg-3
Species	Human
Gene ID (Human)	2078
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human ERG (NP_891548.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P11308
Protein Weight	54kDa
Shipping	Ice bag