

Rabbit Anti-Human Fap Monoclonal Antibody For Fc - Q12884

Item Number: CM0009346



Introduction

-conjugated rabbit anti-human FAP monoclonal antibody for flow cytometry. Specifically detects fibroblast activation protein alpha (FAP α) in human samples. Ideal for cancer, fibrosis, and tissue remodeling research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human FAP mAb
Remarks/Alias	FAP; DPPIV; FAPA; FAPalpha; SIMP; prolyl endopeptidase FAP
Species	Human
GeneID (Human)	2191
GeneID	2191
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-760 of human FAP (NP_004451.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q12884
Protein Weight	27kDa/87kDa
Shipping	Ice bag