

Rabbit Anti-Human Fatp1 Monoclonal Antibody For Intracellular Flow Cytometry - Q6Pcb7

Item Number: CM0012124

Introduction

-conjugated rabbit monoclonal antibody against human FATP1 (SLC27A1), validated for intracellular flow cytometry. Useful for long-chain fatty acid transport and metabolism research. UniProt Q6PCB7.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human FATP1 mAb
Remarks/Alias	FATP; FATP1; ACSVL5; FATP-1
Species	Human
GeneID (human)	376497
GeneID	376497
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 510-646 of human FATP1 (NP_940982.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	Q6PCB7
Protein Weight	52kDa/71kDa
Shipping	Ice bag