


Rabbit Anti-Human FcεR1A Monoclonal Antibody For Flow Cytometry - P12319

Item Number: CM0012614



Introduction

-conjugated rabbit monoclonal antibody targeting human FcεR1α (Fc-epsilon RI-alpha, FCER1A). Validated for flow cytometry (FC) applications. Suitable for allergy, immunology, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human FcεR1α mAb
Remarks/Alias	FCE1A; FcERI; FCERIA
Species	Human
GeneID (Human)	2205
GeneID	2205
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-205 of human FcεR1α (NP_001992.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P12319
Protein Weight	30kDa
Shipping	Ice bag