

Rabbit Anti-Human Fgfr3 Monoclonal Antibody For Flow Cytometry - P22607

Item Number: CM0004927



Introduction

-conjugated rabbit monoclonal antibody against human FGFR3, validated for flow cytometry. FGFR3 is a key regulator of chondrocyte differentiation and skeletal development. Ideal for cellular analysis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human FGFR3 mAb
Remarks/Alias	ACH; CEK2; JTK4; CD333; HSFGR3EX
Species	Human
GeneID (Human)	2261
GeneID	2261
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P22607
Protein Weight	88kDa
Shipping	Ice bag