


Rabbit Anti-Human Flt3/Cd135 Monoclonal Antibody For Flow Cytometry - P36888

Item Number: CM0009648



Introduction

-conjugated rabbit monoclonal antibody targeting human Flt3/CD135. Validated for flow cytometry, this antibody detects the 108-112 kDa receptor tyrosine kinase, essential for hematopoietic progenitor cell studies and leukemia research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Flt3/CD135 mAb
Remarks/Alias	FLT3; CD135; FLK-2; FLK2; STK1; fms related tyrosine kinase 3
Species	Human
Gene ID (Human)	2322
Gene ID	2322
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 27-163 of human Flt3/CD135 (NP_004110.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P36888
Protein Weight	108kDa/112kDa
Shipping	Ice bag