


Rabbit Anti-Human Galectin 3 Mab For Fc (Intra) - P17931

Item Number: CM0008049



Introduction

conjugated rabbit monoclonal antibody specific for human Galectin-3. Ideal for intracellular flow cytometry. Detects a 26 kDa protein. Cross-reacts with human samples. Suitable for immunology, cancer, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Galectin 3 mAb
Remarks/Alias	L31; GAL3; MAC2; CBP35; GALBP; GALIG; LGALS2
Species	Human
Gene ID (Human)	3958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human Galectin 3 (NP_002297.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P17931
Protein Weight	26kDa
Shipping	Ice bag