


## Rabbit Anti-Human Galectin 3 Monoclonal Antibody For Fc (Intra) - P17931

Item Number: CM0012190



### Introduction

-conjugated rabbit monoclonal antibody targeting human Galectin-3 (LGALS3) for intracellular flow cytometry. Recognizes human Galectin-3, ~26kDa. Suitable for immunology, inflammation, and cancer research.

[Learn More](#)

| Parameter        | Value                                        |
|------------------|----------------------------------------------|
| Product Name     | Rabbit anti-Human Galectin 3 mAb             |
| Remarks/Alias    | L31; GAL3; MAC2; CBP35; GALBP; GALIG; LGALS2 |
| Species          | Human                                        |
| GeneID (Human)   | 3958                                         |
| GeneID           | 3958                                         |
| Immunogen        | Recombinant protein of human Galectin 3      |
| Source           | Rabbit                                       |
| Category         | Monoclonal Antibodies                        |
| Application      | FC (intra)                                   |
| Cross Reactivity | Human                                        |
| SWISS            | P17931                                       |
| Protein Weight   | 26kDa                                        |
| Shipping         | Ice bag                                      |