

Rabbit Anti-Human Gata3 Mab For Flow Cytometry - P23771

Item Number: CM0007890



Introduction

-conjugated rabbit monoclonal antibody targeting human GATA3 transcription factor for intracellular flow cytometry. Validated for human reactivity, 48 kDa. Ideal for Th2 cell and macrophage activation studies.

[Learn More](#)

Parameter	Description
Product Name	Rabbit anti-Human GATA3 mAb
Remarks/Alias	HDR; HDRS
Species	Human
GeneID (human)	2625
GeneID	2625
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human GATA3 (NP_001002295.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P23771
Protein Weight	48kDa
Shipping	Ice bag