

Rabbit Anti-Human Glypican 3 (Gpc3) Monoclonal Antibody For Flow Cytometry - P51654

Item Number: CM0009054



Introduction

-conjugated rabbit anti-human Glypican 3 (GPC3) monoclonal antibody, validated for flow cytometry (FC) with human reactivity. GPC3 is a cell surface proteoglycan implicated in Wnt and hedgehog signaling, often investigated in cancer and developmental biology.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Glypican 3 (GPC3) mAb
Remarks/Alias	SGB; DGSX; MXR7; SDYS; SGBS; OCI-5; SGBS1; GTR2-2
Species	Human
Gene ID	2719
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P51654
Protein Weight	66kDa
Shipping	Ice bag