

Rabbit Anti-Human Gp6/Gpvi Monoclonal Antibody For Flow Cytometry - Q9Hcn6

Item Number: CM0008933



Introduction

-conjugated rabbit monoclonal antibody specific for human platelet glycoprotein VI (GP6/GPVI), designed for flow cytometry. Key tool for investigating collagen-induced platelet activation, thrombosis, and hemostasis pathways.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human GP6/GPVI mAb
Remarks/Alias	GPVI; GPVI; BDPLT11
Species	Human
GeneID (Human)	51206
GeneID	51206
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-267 of human GP6/GPVI (NP_057447.5).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Human
SWISS	Q9HCN6
Protein Weight	37kDa
Shipping	Ice bag