


Rabbit Anti-Human Gsk3B Monoclonal Antibody For Flow Cytometry - P49841

Item Number: CM0012809



Introduction

-conjugated rabbit anti-human GSK3 β monoclonal antibody designed for intracellular flow cytometry detection of human GSK3 β , a key kinase in Wnt signaling, glucose metabolism, and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human GSK3 β mAb
Remarks/Alias	GSK3B
Species	Human
GeneID (Human)	2932
GeneID	2932
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 321-420 of human GSK3 β (NP_001139628.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P49841
Protein Weight	47kDa
Shipping	Ice bag