

Rabbit Anti-Human Her2/ErbB2 Monoclonal Antibody For Fc - P04626

Item Number: CM0006120



Introduction

-conjugated rabbit monoclonal antibody targeting human HER2/ErbB2, validated for flow cytometry (FC). Recognizes the extracellular domain (aa 23-652) of this receptor tyrosine kinase, a key oncogene overexpressed in breast and other cancers. Ideal for oncology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human HER2/ErbB2 mAb
Remarks/Alias	NEU; NGL; HER2; TKR1; CD340; HER-2; VSCN2; MLN 19; MLN-19; c-ERB2; c-ERB-2; HER-2/neu; p185(erbB2)
Species	Human
GeneID (Human)	2064
GeneID	2064
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-652 of human HER2/ErbB2 (NP_004439.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P04626
Protein Weight	138kDa
Shipping	Ice bag