


Rabbit Anti-Human Her2/ErbB2 Monoclonal Antibody For Flow Cytometry - P04626

Item Number: CM0012731



Introduction

-conjugated rabbit monoclonal antibody targeting human HER2/ErbB2. Validated for flow cytometry. Immunogen: recombinant extracellular domain (aa 23-652). Use in cancer biology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human HER2/ErbB2 mAb
Remarks/Alias	NEU; NGL; HER2; TKR1; CD340; HER-2; VSCN2; MLN 19; MLN-19; c-ERB2; c-ERB-2; HER-2/neu; p185(erbB2)
Species	Human
Gene ID (Human)	2064
Gene ID	2064
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-652 of human HER2/ErbB2 (NP_004439.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P04626
Protein Weight	138kDa
Shipping	Ice bag