


Rabbit Anti-Human Icam-1/Cd54 Monoclonal Antibody For Flow Cytometry - P05362

Item Number: CM0005769



Introduction

-conjugated rabbit anti-human ICAM-1/CD54 monoclonal antibody (clone not specified) validated for flow cytometry (FC). Detects human ICAM-1, a cell adhesion molecule and rhinovirus receptor, with a predicted molecular weight of 58 kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ICAM-1/CD54 mAb
Remarks/Alias	BB2; CD54; P3.58
Species	Human
Gene ID (Human)	3383
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-480 of human ICAM-1/CD54 (NP_000192.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P05362
Protein Weight	58kDa
Shipping	Ice bag