


Rabbit Anti-Human Icam2/Cd102 Mab For Fc - P13598

Item Number: CM0005215



Introduction

-conjugated rabbit anti-human ICAM2/CD102 monoclonal antibody for flow cytometry. Recombinant immunogen (aa 25-223). Validated for human samples. Ideal for immune cell adhesion and migration studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ICAM2/CD102 mAb
Remarks/Alias	CD102
Species	Human
GeneID (Human)	3384
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-223 of human ICAM2/CD102 (NP_000864.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13598
Protein Weight	31kDa
Shipping	Ice bag