


Rabbit Anti-Human Icam2/Cd102 Monoclonal Antibody For Fc - P13598

Item Number: CM0007992



Introduction

-conjugated rabbit monoclonal antibody against human ICAM2/CD102, validated for flow cytometry. Suitable for studying leukocyte adhesion, immune surveillance, and cell migration. Specific for human ICAM2.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human ICAM2/CD102 mAb
Remarks/Alias	ICAM2; CD102; intercellular adhesion molecule 2
Species	Human
Gene ID (Human)	3384
Gene ID	3384
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-223 of human ICAM2/CD102 (NP_000864.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P13598
Protein Weight	30kDa
Shipping	Ice bag