


Rabbit Anti-Human Icam2/Cd102 Monoclonal Antibody For Flow Cytometry - P13598

Item Number: CM0004209



Introduction

-conjugated rabbit monoclonal antibody targeting human ICAM2/CD102 (SWISS: P13598). Validated for flow cytometry, suitable for leukocyte adhesion and endothelial cell migration research.

[Learn More](#)

-	-
Product Name	Rabbit anti-Human ICAM2/CD102 mAb
Remarks/Alias	CD102
Species	Human
GeneID (human)	3384
GeneID	3384
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-223 of human ICAM2/CD102 (NP_000864.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Human
SWISS	P13598
Protein Weight	31kDa
Shipping	Ice bag