


Rabbit Anti-Human Ido1 Monoclonal Antibody For Flow Cytometry - P14902

Item Number: CM0012343



Introduction

-conjugated rabbit monoclonal antibody against human IDO1 (Indoleamine 2,3-dioxygenase 1) for intracellular flow cytometry. Specific for human samples, 45 kDa, UniProt P14902.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IDO1 mAb
Remarks/Alias	IDO; INDO; IDO-1
Species	Human
GeneID	3620
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-403 of human IDO1 (NP_002155.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P14902
Protein Weight	45kDa
Shipping	Ice bag