

Rabbit Anti-Human Ifn-Alpha/Beta R2 Monoclonal Antibody For Flow Cytometry - P48551

Item Number: CM0009622



Introduction

-conjugated rabbit monoclonal antibody targeting human IFN-alpha/beta R2 (IFNAR2) for flow cytometry. Ideal for studying type I interferon signaling, JAK-STAT pathway, and immune response. Cross-reactive with human.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IFN-alpha/beta R2 mAb
Remarks/Alias	IFN-R; IMD45; IFNABR; IFNARB; IFN-R-2; IFN-alpha-REC; IFNAR2
Species	Human
GeneID (Human)	3455
GeneID	3455
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P48551
Protein Weight	58kDa/37kDa/27kDa
Shipping	Ice bag