

Rabbit Anti-Human Ifn-Gamma Mab For Fc (Intra) - P01579

Item Number: CM0008713



Introduction

High-quality rabbit monoclonal antibody against human interferon gamma (IFN- γ), conjugated to , validated for intracellular flow cytometry (FC intra). Ideal for immune response and cancer immunotherapy research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IFN-gamma mAb
Remarks/Alias	IFG; IFI; IMD69
Species	Human
GeneID (Human)	3458
GeneID	3458
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 24-166 of human IFN-gamma (NP_000610.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P01579
Protein Weight	19 kDa
Shipping	Ice bag