

Rabbit Anti-Human Ifn-Gamma Monoclonal Antibody For Flow Cytometry (Intracellular) - P01579

Item Number: CM0010770



Introduction

-conjugated rabbit monoclonal antibody targeting human IFN-gamma, validated for intracellular flow cytometry. Suitable for studying Th1 immune responses, cytokine signaling, and antitumor immunity.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IFN-gamma mAb
Remarks/Alias	IFG; IFI; IMD69
Species	Human
GeneID	3458 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-166 of human IFN-gamma (NP_000610.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P01579
Protein Weight	19kDa
Shipping	Ice bag