


Rabbit Anti-Human Ig K Light Chain Mab For Fc - P01834

Item Number: CM0005680



Introduction

-conjugated rabbit anti-human Ig κ light chain monoclonal antibody for flow cytometry (FC). Targets the constant region (P01834, ~12 kDa). Reliable tool for B-cell immunology, antibody profiling, and humoral immunity research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Ig κ Light Chain mAb
Remarks/Alias	Km; HCAK1; IGKCD
Species	Human
Gene ID (Human)	3514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-107 of human Ig κ Light Chain (P01834)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01834
Protein Weight	12 kDa
Shipping	Ice bag