


Rabbit Anti-Human Ig K Light Chain Mab - P01834

Item Number: CM0004946



Introduction

-conjugated rabbit monoclonal antibody targeting human Ig κ light chain constant region (IGKC), validated for flow cytometry. Suitable for B cell and immunoglobulin research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Ig κ Light Chain mAb
Remarks/Alias	IGKC, HCAK1, IGKCD, Km
Species	Human
Gene ID (Human)	3514
Gene ID	3514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-107 of human Ig κ Light Chain (P01834)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01834
Protein Weight	11KDa
Shipping	Ice bag