


Rabbit Anti-Human Ig K Light Chain Monoclonal Antibody For Fc - P01834

Item Number: CM0009302



Introduction

-conjugated rabbit monoclonal antibody against human Ig κ light chain constant region (IGKC) for flow cytometry. Specific for human, ~11 kDa, validated in FC. IVD raw material for research use.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Ig κ Light Chain mAb
Remarks/Alias	IGKC, HCAK1, IGKCD, Km
Species	Human
GeneID (human)	3514
GeneID	3514
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-107 of human Ig κ Light Chain (P01834)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01834
Protein Weight	11KDa
Shipping	Ice bag