


Rabbit Anti-Human Ig K Light Chain Monoclonal Antibody For Flow Cytometry - P01834

Item Number: CM0012602



Introduction

-conjugated rabbit monoclonal antibody targeting human immunoglobulin kappa light chain constant region (P01834). Validated for flow cytometry (FC) applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human Ig κ Light Chain mAb
Remarks/Alias	Km; HCAK1; IGKCD
Species	Human
GeneID (human)	3514
GeneID	3514
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-107 of human Ig κ Light Chain(P01834).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01834
Protein Weight	12 kDa
Shipping	Ice bag