


Rabbit Anti-Human IgA1 Monoclonal Antibody For Fc - P01876

Item Number: CM0011025



Introduction

-conjugated rabbit monoclonal antibody against human IgA1 constant region. Validated for flow cytometry (FC). Recognizes human IgA1. Suitable for detection of IgA1-positive cells in mucosal immunity research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IgA1 mAb
Remarks/Alias	IgA1
Species	Human
GeneID (human)	3493
GeneID	3493
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 99-353 of human IgA1.
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01876
Protein Weight	38kDa/42KDa
Shipping	Ice bag