

Rabbit Anti-Human IgE Monoclonal Antibody For Fc - P01854

Item Number: CM0011731



Introduction

-conjugated rabbit monoclonal antibody targeting human IgE constant region. Validated for flow cytometry (FC), ideal for B cell immunophenotyping and IgE receptor studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IgE mAb
Remarks/Alias	IgE
Species	Human
Gene ID (Human)	3497
Gene ID	3497
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 104-428 of human IgE (P01854).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01854
Protein Weight	47kDa / 54kDa / 59kDa
Shipping	Ice bag