

Rabbit Anti-Human IgE Monoclonal Antibody For Fc - P01854

Item Number: CM0010896



Introduction

-conjugated rabbit monoclonal antibody against human IgE heavy chain, validated for flow cytometry. Suitable for investigating IgE-mediated hypersensitivity, allergy, and immune cell activation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IgE mAb
Remarks/Alias	IgE
Species	Human
GeneID (Human)	3497
GeneID	3497
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 104-428 of human IgE (P01854)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01854
Protein Weight	47kDa/54kDa/59kDa
Shipping	Ice bag