


Rabbit Anti-Human IgG (Fc) Monoclonal Antibody For Fc - P01857

Item Number: CM0006203



Introduction

-conjugated rabbit monoclonal antibody targeting human IgG (Fc) region. Validated for flow cytometry (FC). Ideal for IgG detection and Fc receptor studies. SWISS: P01857.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IgG (Fc) mAb
Species	Human
Gene ID	3500 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-330 of human IgG1 (Fc) (P01857)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01857
Protein Weight	36 kDa
Shipping	Ice bag