


Rabbit Anti-Human IgG (Fc) Monoclonal Antibody For Flow Cytometry - P01857

Item Number: CM0012229



Introduction

-conjugated rabbit monoclonal antibody targeting human IgG Fc region (SWISS: P01857). Validated for flow cytometry; specific for human IgG1 Fc. Ideal for immune cell profiling, antibody detection, and functional studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IgG (Fc) mAb
Species	Human
Gene ID	3500
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 100-330 of human IgG1 (Fc) (P01857).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01857
Protein Weight	36kDa
Shipping	Ice bag