


Rabbit Anti-Human IgG (Fc) Monoclonal Antibody For Flow Cytometry - P01857

Item Number: CM0007764



Introduction

-conjugated rabbit monoclonal antibody targeting human IgG (Fc) region, validated for flow cytometry (FC). Useful for detecting membrane IgG or immune complexes in immunological research. P01857.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IgG (Fc) mAb
Species	Human
GeneID (Human)	3500
GeneID	3500
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-330 of human IgG (Fc) (P01857).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01857
Shipping	Ice bag