

Rabbit Anti-Human IL-10Ra/Cd210 Mab For Fc - Q13651

Item Number: CM0009171



Introduction

-conjugated rabbit monoclonal antibody targeting human IL-10RA/CD210 (IL-10 receptor alpha) specifically validated for flow cytometry. Ideal for studying anti-inflammatory signaling and immune cell phenotyping.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL-10RA/CD210 mAb
Remarks/Alias	CD210; IL10R; CD210a; CDW210A; HIL-10R; IL-10R1
Species	Human
GeneID (human)	3587
GeneID	3587
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-235 of human IL-10RA/CD210 (NP_001549.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13651
Protein Weight	63 kDa
Shipping	Ice bag