


Rabbit Anti-Human IL-17Rb Monoclonal Antibody For Flow Cytometry - Q9Nrm6

Item Number: CM0012647



Introduction

Rabbit monoclonal antibody against human IL-17RB, conjugated to ABflo 488. Validated for flow cytometry (FC). Targets the IL-17B/E receptor involved in hematopoietic cell growth and differentiation. Research use only.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human IL-17RB mAb
Remarks/Alias	IL17RB; CRL4; EVI27; IL17BR; IL17RH1; interleukin-17 receptor B
Species	Human
GeneID	55540 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-289 of human IL-17RB (NP_061195.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9NRM6
Protein Weight	31kDa/55kDa
Shipping	Ice bag