


Rabbit Anti-Human IL-2 Monoclonal Antibody For Intracellular Flow Cytometry - P60568

Item Number: CM0011709



Introduction

Rabbit monoclonal antibody against human IL-2 (TCGF), conjugated to a fluorophore for intracellular flow cytometry. Suitable for studying immune response, T-cell growth, and tolerance. Reacts with human.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL-2 mAb
Remarks/Alias	IL-2; TCGF; lymphokine
Species	Human
Gene ID (Human)	3558
Gene ID	3558
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-153 of human IL-2 (NP_000577.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P60568
Protein Weight	18kDa
Shipping	Ice bag