


Rabbit Anti-Human IL-4 Monoclonal Antibody For Intracellular Flow Cytometry - P05112

Item Number: CM0010600



Introduction

-conjugated rabbit monoclonal antibody targeting human interleukin-4 (IL-4). Recombinant immunogen; suitable for intracellular flow cytometry on human samples. Ideal for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL-4 mAb
Remarks/Alias	BSF1; IL-4; BCGF1; BSF-1; BCGF-1
Species	Human
Gene ID	3565 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-153 of human IL-4 (NP_000580.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P05112
Protein Weight	17 kDa
Shipping	Ice bag