


Rabbit Anti-Human IL-5 Monoclonal Antibody For Flow Cytometry (Intracellular) - P05113

Item Number: CM0010232



Introduction

-conjugated rabbit monoclonal antibody targeting human interleukin-5 (IL-5). Validated for intracellular flow cytometry. Suitable for research on eosinophil biology, allergic inflammation, and immune responses.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL-5 mAb
Remarks/Alias	EDF; IL-5; TRF; interleukin-5
Species	Human
Gene ID (Human)	3567
Gene ID	3567
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P05113
Protein Weight	15kDa
Shipping	Ice bag