


Rabbit Anti-Human IL-5 Monoclonal Antibody For Flow Cytometry - P05113

Item Number: CM0004687



Introduction

-conjugated rabbit monoclonal antibody targeting human IL-5, a key cytokine for eosinophil regulation. Validated for intracellular flow cytometry (FC) in immunology studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL-5 mAb
Remarks/Alias	EDF; IL-5; TRF; interleukin-5
Species	Human
GeneID	3567
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-134 of human IL-5 (NP_000870.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P05113
Protein Weight	15kDa
Shipping	Ice bag