


Rabbit Anti-Human IL-6R/Cd126 Monoclonal Antibody For Flow Cytometry (Intracellular) - P08887

Item Number: CM0009144



Introduction

-conjugated rabbit anti-human IL-6R/CD126 monoclonal antibody for intracellular flow cytometry. Recombinant extracellular domain immunogen; cross-reacts with human IL-6R. Useful for investigating IL-6 classic and trans-signaling pathways in immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL-6R/CD126 mAb
Remarks/Alias	IL6Q; gp80; CD126; HIES5; IL-6R; IL6RA; IL6RQ; IL-1Ra; IL-6RA; IL6QTL; IL-6R-1
Species	Human
GeneID (Human)	3570
GeneID	3570
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-358 of human IL-6R/CD126 (NP_000556.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P08887
Protein Weight	52kDa
Shipping	Ice bag