


Rabbit Anti-Human IL18R1 Monoclonal Antibody For Flow Cytometry - Q13478

Item Number: CM0005692



Introduction

-conjugated rabbit monoclonal antibody against human IL18R1 (CD218a), validated for flow cytometry. Ideal for investigating IL18-mediated pathways, Th1 cell responses, and inflammatory disease research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL18R1/IL-18 R alpha mAb
Remarks/Alias	CD218a; IL18RA; IL1RRP; CDw218a; IL-1Rrp; IL-18Ralpha; IL18Ralpha2; IL-18R-alpha
Species	Human
GeneID (Human)	8809
GeneID	8809
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-329 of human IL18R1/IL-18 R alpha (NP_003846.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13478
Protein Weight	62kDa
Shipping	Ice bag