


Rabbit Anti-Human IL23R Monoclonal Antibody For Flow Cytometry - Q5Vwk5

Item Number: CM0008469



Introduction

-conjugated rabbit monoclonal antibody against human IL23R (PSORS7). Validated for flow cytometry (FC). Ideal for studying IL-23 signaling in immune cells, autoimmune diseases, and tumor immunology.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL23R mAb
Remarks/Alias	PSORS7
Species	Human
GeneID (Human)	149233
GeneID	149233
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 124-353 of human IL23R (NP_653302.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q5VWK5
Protein Weight	41kDa/43kDa/72kDa
Shipping	Ice bag