


Rabbit Anti-Human IL8 Monoclonal Antibody For Flow Cytometry (Intracellular) - P10145

Item Number: CM0012469



Introduction

-conjugated rabbit monoclonal antibody against human IL8 for intracellular flow cytometry. Validated for human samples, targeting the chemokine IL8 involved in inflammatory responses. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL8 mAb
Remarks/Alias	IL8; NAF; GCP1; LECT; LUCT; NAP1; GCP-1; LYNAP; MDNCF; MONAP; NAP-1; SCYB8
Species	Human
Gene ID (Human)	3576
Gene ID	3576
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-99 of human IL8 (NP_000575.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P10145
Protein Weight	11kDa
Shipping	Ice bag