


Rabbit Anti-Human IL8 Monoclonal Antibody For Flow Cytometry - P10145

Item Number: CM0006079



Introduction

-conjugated rabbit monoclonal antibody targeting human IL8 (Interleukin-8), validated for intracellular flow cytometry. Ideal for studying chemokine-driven inflammation, with a protein weight of 11kDa.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL8 mAb
Remarks/Alias	IL8; NAF; GCP1; LECT; LUCT; NAP1; GCP-1; LYNAP; MDNCF; MONAP; NAP-1; SCYB8
Species	Human
GeneID	3576
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P10145
Protein Weight	11kDa
Shipping	Ice bag