


Rabbit Anti-Human IL8 Monoclonal Antibody For Intracellular Flow Cytometry - P10145

Item Number: CM0004351



Introduction

Rabbit anti-Human IL8 monoclonal antibody validated for intracellular flow cytometry. Detects human IL8, a key chemokine in inflammation and immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human IL8 mAb
Remarks/Alias	IL8; NAF; GCP1; LECT; LUCT; NAP1; GCP-1; LYNAP; MDNCF; MONAP; NAP-1; SCYB8
Species	Human
Gene ID (Human)	3576
Gene ID	3576
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P10145
Protein Weight	11kDa
Shipping	Ice bag